

Electrical Compatible Products

Below is a list of selected key Eve Taylor® products compatible with various devices to achieve outstanding treatment room results.

Microdermabrasion	Galvanic	Microcurrent & Faradic	Ultrasonic	High Frequency	LED Light Therapy
<p>Prepare - Clear Cleanse Purifying Wash</p> <p>Hydrate & treat - Hydrating Serum Soothing Masque Bio Peel-off Masque Ultra Soothing Toner</p> <p>Nourish & protect - Refining Eye Gel Ultra Soothing Moisture Cream Moisturising Solar Shield SPF25 Daily Skin Defence SPF50 Seal & Protect Moisturising Lip Balm SPF10</p>	<p>Desincrustation - Desincrustation Gel (-)</p> <p>Iontophoresis - Pigmented skin: Pigment Response Serum (+) Brightening Serum (+)</p> <p>Oily/ breakout skin: Sebum Response Serum (+) Balancing Masque (+) Purifying Masque (+)</p> <p>Mature/ageing Skin: Oxygen Response Serum (+) Firming Serum (+) Timeline Response Serum (+) Anti-oxidant Masque (+)</p> <p>Dry/dehydrated skin: Hydrating Serum (+) Oxygen Response Serum (+)</p> <p>Rosacea/reddened skin: Soothing Masque mixed with Colloidal Activator (+)</p> <p>All skins: Electro Gel (+/-) Soothing Masque (+) Anti-oxidant Masque (+)</p>	<p>Microcurrent - Lifting: Firming Serum Electro Gel</p> <p>Specific iontophoresis/ feathering: Timeline Response Serum Hydrating Serum Refining Eye Gel</p> <p>Faradic Firming Serum Electro Gel</p>	<p>Ultrasonic exfoliation - Dynamic Resurfacing Cleanser Micro-fine Daily Exfoliant Ultra Soothing Toner</p> <p>Ultrasonic clarifying / decongesting - Desincrustation Gel (-) Electro Gel (+/-) with Hydrolat</p> <p>Ultrasonic ionisation - Pigmented skin: Pigment Response Serum (+) Brightening Serum (+)</p> <p>Oily/breakout skin: Sebum Response Serum (+) Balancing Masque (+) Purifying Masque (+)</p> <p>Mature/ageing Skin: Firming Serum (+) Anti-oxidant Masque (+)</p> <p>Dry/dehydrated skin: Hydrating Serum (+) Oxygen Response Serum (+) Soothing Masque (+)</p> <p>Rosacea/reddened skin: Soothing Masque mixed with Colloidal Activator (+)</p>	<p>Indirect high frequency/ vienesse massage - Organic Massage Cream Revitalising Massage Oil</p> <p>Direct high frequency - Sebum Response Serum, (with dry gauze over top)</p> <p>Purifying Spot Gel (applied directly to breakouts prior to 'sparking' with fulgerator)</p>	<p>Pigmented skin: Pigment Response Serum Brightening Serum Anti-oxidant Masque</p> <p>Oily/breakout skin: Sebum Response Serum Purifying Spot Gel</p> <p>Mature/ageing Skin: Firming Serum Rejuvenating Aromatic Serum Timeline Response Serum Anti-oxidant Masque</p> <p>Dry/dehydrated skin: Hydrating Serum Nourishing Serum Oxygen Response Serum Anti-oxidant Masque</p>